




## Product Specification

Urine container, 70 ml, (LxØ): 55 x 44 mm,  
graduated, PP, transparent, with paper label

Page 1



### Product description

Order number	75.9922.745
Product description	Urine container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: green, 250 piece(s)/bag

### Product characteristics

Label/ Print	with paper label
Colour of print/label	green
Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	70 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

Urine container, 70 ml, (LxØ): 55 x 44 mm,  
graduated, PP, transparent, with paper label

Page 2

### Size

Sample volume	70 ml
Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm
Length of product	55 mm
Length excluding cap	55 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	8000
Depth of case	548 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0834 cbm
Weight of product	0.0115 kg
Weight of case	6.6 kg
EAN of inner box	4038917152761
EAN case	4038917027809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**