



Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag

Page 1



Product description

| | |
|---------------------|--|
| Order number | 10.252.028 |
| Product description | Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15.3 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, with plastic label, sterile, 100 piece(s)/bag |

Product characteristics

| | |
|-----------------------|--|
| Type of collection | Aspiration |
| Type of preparation | without preparation |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/yellow |
| Graduation | yes |
| Cap | Luer screw cap |
| Closure type | screw cap |
| Light protection | no |
| Design | Luer lock |
| Application | Urine tube |
| Example of use | Urine collection, transport, and storage |
| Base shape | round base |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag

Page 2

| | |
|----------------------------|-----------|
| Volume of work | 10 ml |
| Recommended centrifugation | 400 x g |
| Centrifugation time | 5 Min |
| Connection | Luer lock |

Size

| | |
|----------------------------------|---------|
| Sample volume | 10 ml |
| Diameter | 15.3 mm |
| Length including cap | 117 mm |
| Length excluding cap | 102 mm |
| Length including cap and plunger | 127 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|------------|
| Minimum order qty. | 100 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 100 |
| Piece(s) / pallet | 7000 |
| Case volume | 0.0229 cbm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15.3 mm, 100 piece(s)/bag

| | |
|-------------------|---------------|
| Weight of product | 0.0149 kg |
| Weight of case | 1.54 kg |
| EAN of inner box | 4038917104494 |
| EAN case | 4038917075671 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

