



Product description

| | |
|---------------------|--|
| Order number | 11.1240.100 |
| Product description | Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag |

Product characteristics

| | |
|--------------------|-------------------------------|
| Type of collection | Vacuum |
| Light protection | no |
| Design | With wings on the tube holder |
| Application | Urine transfer |
| Suitable for | V-Monovette® Urine |
| Connection | Vacuum tube |

Size

| | |
|----------------|-------|
| Length of tube | 10 cm |
|----------------|-------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product yellow

Purity & certification

| | |
|-------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE supplier |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 30000 |
| Depth of box | 300 mm |
| Width of box | 250 mm |
| Height of box | 10 mm |
| Depth of case | 380 mm |
| Width of case | 350 mm |
| Height of case | 170 mm |
| Case volume | 0.0226 cbm |
| Weight of product | 0.002 kg |
| Weight of case | 1.8 kg |
| EAN of inner box | 6947462401157 |
| EAN case | 6947462401423 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

