



Product description

Order number	93.875.281
Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, red, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 95 mm, format: 9 x 9, for 81 collection tubes, for CryoPure tubes 3.5?5.0 ml internal and external thread, 5 piece(s)/bag

Product characteristics

Graduation	no
Cap	slip-on lid with ventilation function
Closure type	slip-on lid
Design	with numerical coding at each aperture
Suitable for	CryoPure tubes 3.5?5.0 ml internal and external thread
Storage capacity	81
Format	9 x 9

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Width of product	132 mm
Height of product	95 mm
Length of product	132 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	red
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	400
Depth of case	305 mm
Width of case	305 mm
Height of case	508 mm
Case volume	0.0472 cbm
Weight of product	0.2437 kg
Weight of case	5 kg
EAN of inner box	4038917122948
EAN case	4038917076647

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

