



Product description

Order number	83.3941.502
Product description	Filtropur V50, Vacuum filtration unit, 500 ml, screw cap, membrane: PES, Ø membraner: 75 mm, pore size: 0.1 µm, for removing mycoplasmas, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Example of use	for removing mycoplasmas
Filter	yes
Surface area of filter	38.04 cm ²
Pore size	0,1
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Filtropur V50, Vacuum filtration unit, 500 ml, PES,
0.1 µm

Page 2

Size

Sample volume	500 ml
Filtration volume	500 ml
Height of product	24.6 cm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Membrane material	Polyethersulfone (PES)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	12
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	288
Depth of case	508 mm
Width of case	391 mm
Height of case	285 mm
Case volume	0.0566 cbm
Weight of product	0.1724 kg
Weight of case	3.8 kg
EAN of inner box	4038917969628
EAN case	4038917969635

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

