



Product description

Order number	83.3942.501
Product description	Filtropur V100, Vacuum filtration unit, 1,000 ml, screw cap, membrane: PES, Ø membraner: 91 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Example of use	for sterile filtration
Filter	yes
Surface area of filter	58.93 cm ²
Pore size	0,22
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Filtropur V100, Vacuum filtration unit, 1,000 ml,
PES, 0.22 µm**

Page 2

Size

Sample volume	1.000.00 ml
Filtration volume	1000 ml
Height of product	31 cm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Membrane material	Polyethersulfone (PES)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	12
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	12
Piece(s) / pallet	192
Depth of case	545 mm
Width of case	425 mm
Height of case	350 mm
Case volume	0.0811 cbm
Weight of product	0.26 kg
Weight of case	5.42 kg
EAN of inner box	4038917969123
EAN case	4038917969130

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-20

This document was prepared by EDP
support and is valid without signature.



SARSTEDT