



**Product Specification**  
**lumox® multiwell, Cell culture plate, with foil**  
**base, 24 well, 4 piece(s)**

**Page 1**



**Product description**

Order number	94.6000.014
Product description	lumox® multiwell, Cell culture plate, with foil base, 24 well, frame: black, frame: PS, for Adherent cells, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Closure type	slip-on lid
Design	with foil base
Cell type	Adherent cells
Base shape	flat base
Number of wells	24
Growth area	1.9 cm <sup>2</sup>
Volume of work	500 µl - 1.500 µl
Vessel type	Plate
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**lumox® multiwell, Cell culture plate, with foil  
base, 24 well, 4 piece(s)**

**Page 2**

**Size**

Width of product	85.48 mm
Height of product	23.4 mm
Length of product	127.76 mm

**Material & colours**

Product material	lumox®
Frame material	Polystyrene (PS)
Colour of frame	black

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	4
Type of smallest subpackaging	blister
Piece(s) / inner box	4
Piece (s) / outer case	16
Piece(s) / pallet	1024
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0653 kg
Weight of case	1.93 kg
EAN of inner box	4038917254571
EAN case	4038917254564

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

