



Product Specification
**Centrifuge, SMC 6 V, supply voltage: 115 V,
swing-out rotor 6-fold**

Page 1

Product description

| | |
|---------------------|--|
| Order number | 90.184.401 |
| Product description | Centrifuge, SMC 6 V, supply voltage: 115 V, incl. swing-out rotor 6-fold and tube holder 1-fold for tubes up to 100 mm length, for separation of serum/plasma from the blood cells, swing-out rotor 6-fold |

Product characteristics

| | |
|---------------------|---|
| Design | swing-out rotor 6-fold |
| Suitable for | separation of serum/plasma from the blood cells |
| Number of positions | 6 |
| Supply voltage | 115 V |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 36 |
| Weight of product | 5.7995 kg |
| Weight of case | 5.8 kg |
| EAN of inner box | 4038917534420 |
| EAN case | 4038917534413 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

