



Product description

Order number	92.2620.301
Product description	length: 1.5 m, material: Bioprene, transparent, for IVARO, 1 piece(s)/bag

Product characteristics

Suitable for	IVARO
--------------	-------

Size

Diameter	2.4 mm
Length of product	1.5 m

Material & colours

Product material	Bioprene
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	10
Depth of case	280 mm
Width of case	220 mm
Height of case	100 mm
Case volume	0.0062 cbm
Weight of product	0.0005 kg
EAN of inner box	4038917584234
EAN case	4038917584227

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

