



Product Specification
**Weigh tray, 250 ml, (LxW): 128 x 128 mm, PVC,
white**

Page 1



Product description

Order number	71.9923.210
Product description	Weigh tray, 250 ml, (LxW): 128 x 128 mm, material: PVC, white, 100 piece(s)/bag

Size

Sample volume	250 ml
Max. volume	300 ml
Width of product	128 mm
Length of product	128 mm

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Weigh tray, 250 ml, (LxW): 128 x 128 mm, PVC,
white**

Page 2

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	28000
Depth of case	458 mm
Width of case	373 mm
Height of case	238 mm
Case volume	0.0407 cbm
Weight of product	0.005 kg
Weight of case	5.6 kg
EAN of inner box	4038917153645
EAN case	4038917038225

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT