



## Product description

Order number	84.2620.003
Product description	Adhesive label, (LxW): 30 x 20 mm, 3,500 piece(s)/roll

## Size

Width of product	20 mm
Length of product	30 mm

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order qty.	7000
Type of smallest subpackaging	roll
Piece(s) / inner box	3500
Piece (s) / outer case	7000
Piece(s) / pallet	2268000
Depth of case	503 mm
Width of case	375 mm
Height of case	155 mm
Case volume	0.0292 cbm
Weight of product	0.0001 kg
Weight of case	0.73 kg
EAN of inner box	4038917900362
EAN case	4038917900379

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

