



### **Product description**

Order number	92.181.611
Product description	Rotor, for 22 collection tubes

### **Product characteristics**

Storage capacity	22
------------------	----

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	40
Case volume	0.0052 cbm
Weight of product	1.3595 kg
Weight of case	1.36 kg
EAN of inner box	4038917156257
EAN case	4038917156240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

