



## **Product description**

Order number	97.8710.550
Product description	Module linfusion heater, for warming infusion solutions, tubes, instruments and contrast agent in the SAHARA-III

## **Product characteristics**

Application	Warming
-------------	---------

### **Size**

Width of product	4.5 cm
Length of product	6 cm

## **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0039 cbm
Weight of product	0.03 kg
Weight of case	0.22 kg
EAN of inner box	4038917968379
EAN case	4038917968386

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

