



Product description

Order number	76.1340.500
Product description	Lid, for 30-ml medicine cup, material: PE, white, 750 piece(s)/bag

Product characteristics

Example of use	Oral administration of medications in liquid or tablet form
Suitable for	30-ml medicine cup

Material & colours

Product material	Polyethylene (PE)
Colour of product	white
Colour of cap	white

Purity & certification

Product category	Medical device
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



