



Product Specification
**Micro sample tube Lithium heparin LH, 1.3 ml,
push cap, EU**

Page 1



Product description

| | |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 41.1503.015 |
| Product description | Micro sample tube Lithium heparin LH, venous blood collection, preparation: Lithium heparin, nominal volume: 1.3 ml, Ø: 10.8 mm, push cap, cap: orange, colour code EU, with paper label, label/print: orange, skirted conical base, 100 piece(s)/bag |

Product characteristics

| | |
|---------------------------|-------------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin |
| Preparation concentration | 25 IU/ml blood |
| Colour code | EU |
| Label/ Print | with paper label |
| Colour of print/label | orange |
| Graduation | yes |
| Cap | flat cap |
| Closure type | push cap |
| Application | venous blood collection |
| Base shape | skirted conical base |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Micro sample tube Lithium heparin LH, 1.3 ml,
push cap, EU**

Page 2

Volume of work 1,3 ml
Recommended centrifugation 2000 x g
Centrifugation time 5 Min

Size

Sample volume 1.3 ml
Diameter 10.8 mm
Length of product 44 mm
Length excluding cap 40 mm

Material & colours

Colour of product transparent
Closure material Low Density Polyethylene (LD-PE)
Colour of cap orange
Tube material Polypropylene (PP)

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0016 kg
Weight of case 1.89 kg
EAN of inner box 4038917310192
EAN case 4038917310185

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

