



**Product description**

Order number 71.9923.211  
Product description Weigh tray, 5 ml, (LxW): 35 x 35 mm, material: PVC, white, 200 piece(s)/bag

**Size**

Sample volume 5 ml  
Max. volume 8 ml  
Width of product 35 mm  
Length of product 35 mm

**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product white

**Purity & certification**

Product category no medical device | IVD  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	308000
Depth of case	271 mm
Width of case	158 mm
Height of case	224 mm
Case volume	0.0096 cbm
Weight of product	0.0007 kg
Weight of case	1.6 kg
EAN of inner box	4038917153652
EAN case	4038917038232

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

