



## Product description

Order number	20.1330
Product description	Microvette® APT 500 EDTA K2E, capillary blood collection, preparation: K2 EDTA, K2E, nominal volume: 500 µl, Ø: 13 mm, with plastic label, cap red, cap, with pierceable membrane, round base, inner receptacle: cylindrical, non-sterile, 100 piece(s)/case

## Product characteristics

Type of collection	capillary
Type of preparation	K2 EDTA
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/red
Cap	membrane cap
Closure type	cap
Application	capillary blood collection
Base shape	round base
Volume of work	250 µl - 500 µl
Shape of inner tube	cylindrical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Sample volume	500 µl
Diameter	13 mm
Length excluding cap	75 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polyethylene (PE)
Colour of cap	red

## Purity & certification

Satisfies the requirement	P 650, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

## Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	54000
Depth of box	166 mm
Width of box	176 mm
Height of box	68 mm
Depth of case	358 mm
Width of case	197 mm
Height of case	186 mm
Case volume	0.0131 cbm
Weight of product	0.006 kg
Weight of case	3.17 kg
EAN of inner box	4038917524735
EAN case	4038917524728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

