



**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 51.9923.820   |
| Product description | Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, skirted conical base, transparent, cap assembled, with paper label, preparation: Boric acid, non-sterile, 500 piece(s)/case |

**Product characteristics**

|                          |                      |
|--------------------------|----------------------|
| Type of preparation      | Boric acid           |
| Label/ Print             | with paper label     |
| Colour of print/label    | black                |
| Closure type             | screw cap            |
| Base shape               | skirted conical base |
| Volume of work           | 25 ml                |
| Centrifugation max (RCF) | 4500 x g             |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile**

**Page 2**

**Size**

|                   |       |
|-------------------|-------|
| Sample volume     | 25 ml |
| Diameter          | 25 mm |
| Length of product | 90 mm |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Colour of cap     | red                |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | non-sterile                   |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 10000         |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 15000         |
| Case volume                   | 0.059 cbm     |
| Weight of product             | 0.009 kg      |
| Weight of case                | 5.18 kg       |
| EAN of inner box              | 4038917154550 |
| EAN case                      | 4038917030090 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**