



**Product description**

Order number	74.5684.001
Product description	Urine measuring system exchange bag, 2 l, sterile, 1 piece(s)/blister

**Product characteristics**

Colour of print/label	blue
Application	Urine drainage
Volume of work	2 l

**Size**

Sample volume	2 l
Max. volume	2 l

**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	4050
Depth of case	395 mm
Width of case	195 mm
Height of case	200 mm
Case volume	0.0154 cbm
Weight of product	0.0451 kg
Weight of case	3.02 kg
EAN of inner box	6955824024432
EAN case	6955824024050

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

