




**Product Specification**  
**Serological pipette, with tip, plugged, 5 ml, sterile,**  
**25 piece(s)/bag**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 86.1253.025   |
| Product description | Serological pipette, with tip, plugged, 5 ml, coding colour: blue, material: PS, sterile, 25 piece(s)/bag |

**Product characteristics**

|                       |                   |
|-----------------------|-------------------|
| Colour of print/label | blue              |
| Graduation            | yes               |
| Design                | with tip, plugged |
| Volume of work        | 5 ml              |
| Coding colour         | blue              |

**Size**

|                   |        |
|-------------------|--------|
| Diameter          | 8 mm   |
| Length of product | 340 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | sterile                 |
| Sterilisation    | Gamma irradiation       |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 25            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 40500         |
| Depth of case                 | 413 mm        |
| Width of case                 | 250 mm        |
| Height of case                | 195 mm        |
| Case volume                   | 0.0201 cbm    |
| Weight of product             | 0.0067 kg     |
| Weight of case                | 4.03 kg       |
| EAN of inner box              | 4038917127257 |
| EAN case                      | 4038917021111 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

