



### **Product description**

Order number 95.1710  
Product description Foam film bag, (LxW): 185 x 230 mm, material: PE

### **Size**

Width of product 230 mm  
Length of product 185 mm

### **Material & colours**

Product material Polyethylene (PE)

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 200  
Piece(s) / inner box 200  
Piece (s) / outer case 200  
Piece(s) / pallet 6400  
Depth of case 458 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Foam film bag, (LxW): 185 x 230 mm, PE**

**Page 2**

Width of case	408 mm
Height of case	205 mm
Case volume	0.0383 cbm
Weight of product	0.0018 kg
Weight of case	0.91 kg
EAN of inner box	4038917087490
EAN case	4038917079648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**