



### Product description

Order number	65.809.306
Product description	Push cap, red, suitable for tubes Ø 12 mm, 1,000 piece(s)/bag

### Product characteristics

Closure type	push cap
Suitable for	tubes Ø 12 mm

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	red

### Purity & certification

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Push cap, red, suitable for tubes Ø 12 mm**

**Page 2**

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	400000
Case volume	0.0382 cbm
Weight of product	0.0006 kg
Weight of case	6.82 kg
EAN of inner box	4038917369206
EAN case	4038917369190

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**