



Product description

Order number	92.180.613
Descrição do artigo	Rotor em bloco, para 24 tubos até 15 mm Ø, para SARMIX® M 2000, 1 unid./caixa de papelão

Product characteristics

Indicado para	SARMIX® M 2000
Capacidade de armazenamento	24

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	não

Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	10
Unidades por palete	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rotor em bloco, para 24 tubos até 15 mm Ø, para
SARMIX® M 2000

Page 2

Volume da embalagem externa	0.0122 cbm
Peso do produto	0.1509 kg
Peso da embalagem externa	1.71 kg
EAN da caixa interna	4038917102018
EAN da embalagem externa	4038917075770

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT