



Product Specification
Deep Well MegaBlock® , 0,5 ml, Isento de DNA/DNase/RNase, isento de pirogênios/endotoxinas, PP

Page 1



Product description

Order number	82.1969.002
Descrição do artigo	Deep Well MegaBlock® , 96 poço, 0,5 ml, fundo redondo, Isento de DNA/DNase/RNase, isento de pirogênios/endotoxinas, material: PP, padrão SBS, cavidades redondas, 7 unid./pacote

Product characteristics

Tipo de fechamento	filme
Versão	padrão SBS
Campo de aplicação	coleta e armazenamento de amostras para diagnóstico
Forma do fundo	fundo redondo
Quantidade de poços	96
Volume de trabalho	1 µl - 0,5 ml
Centrifugação máxima (RCF)	2000 x g
Temperatura de centrifugação	22
Tempo de centrifugação	5 Min
Autoclavável (máx. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Deep Well MegaBlock® , 0,5 ml, Isento de DNA/DNase/RNase, isento de pirogênios/endotoxinas, PP

Page 2

Size

Volume nominal	0.5 ml
Diâmetro	7.2 mm
Largura do produto	85.1 mm
Altura do produto	28.5 mm
Comprimento do produto	127.3 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Padrão de pureza	Isento de DNA/DNase/RNase, isento de pirogênios/endotoxinas
Em lotes	sim

Packaging

Quantidade mínima do pedido	56
Tipo de embalagem interior	pacote
Unidades por caixa interna	7
Unidades por caixa	56
Unidades por palete	4480
Volume da embalagem externa	0.0181 cbm
Peso do produto	0.057 kg
Peso da embalagem externa	3.74 kg
EAN da caixa interna	4038917346320
EAN da embalagem externa	4038917346313

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT