




**Product Specification**  
**TUR-BAG, Urine drainage system, 4 l, without**  
**pump ball, sterile**

**Page 1**



**Product description**

Order number	74.5340.501
Product description	TUR-BAG, Urine drainage system, 4 l, without pump ball, sterile, tube length: 100 cm, 1 piece(s)/blister

**Product characteristics**

Colour of print/label	blue
Design	without pump ball
Application	Urine drainage
Volume of work	4 l

**Size**

Sample volume	4 l
Max. volume	4 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**TUR-BAG, Urine drainage system, 4 l, without  
pump ball, sterile**

**Page 2**

**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

**Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	960
Depth of case	395 mm
Width of case	295 mm
Height of case	260 mm
Case volume	0.0303 cbm
Weight of product	0.1645 kg
Weight of case	3.83 kg
EAN of inner box	6955824024401
EAN case	6955824024029

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-02

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**