




Product Specification
**Petri dish, 60 x 15 mm, transparent, with
ventilation cams**

Page 1



Product description

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------|
| Order number | 82.1194.500 |
| Product description | Petri dish, (ØxH): 60 x 15 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag |

Product characteristics

| | |
|---------------------|-------|
| Minimum temperature | 0 °C |
| Maximum temperature | 80 °C |
| Ventilation cams | yes |

Size

| | |
|-------------------|-------|
| Max. volume | 30 ml |
| Diameter | 60 mm |
| Height of product | 15 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Petri dish, 60 x 15 mm, transparent, with
ventilation cams**

Page 2

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|-------------------|-------------------------|
| Product category | no medical device IVD |
| Latex-free | yes |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Safety data sheet | no |
| Batched | yes |

Packaging

| | |
|----------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | tube bag |
| Piece(s) / inner box | 20 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 12000 |
| Depth of case | 317 mm |
| Width of case | 317 mm |
| Height of case | 365 mm |
| Case volume | 0.0367 cbm |
| Weight of product | 0.0074 kg |
| Weight of case | 4.49 kg |
| EAN of inner box | 4038917347891 |
| EAN case | 4038917347884 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-03

This document was prepared by EDP
support and is valid without signature.



SARSTEDT