



Product description

Order number 71.9924.515
Product description Forceps, height: 11 mm, green, PS

Product characteristics

Application Collection of diagnostic specimens

Size

Width of product 28 mm
Height of product 11 mm
Length of product 108 mm

Material & colours

Product material Polystyrene (PS)
Colour of product green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	3000
Piece(s) / inner box	3000
Piece (s) / outer case	3000
Piece(s) / pallet	90000
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of case	12 kg
EAN of inner box	4038917155281
EAN case	4038917398398

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

