



Product description

Order number	87.7709.029
Product description	Tube stripper AST-2, battery-powered, incl. transport case

Product characteristics

Application	Obtaining blood and blood products
Supply voltage	240 V
Power consumption during operation	20 W
Frequency	50 Hz - 60 Hz
Battery type	Nickel metal hydride battery

Size

Width of product	39 mm
Height of product	62 mm
Length of product	225 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	not regulated
CE certified	CE - manufacturer's self-declaration
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	84
Depth of case	400 mm
Width of case	360 mm
Height of case	130 mm
Case volume	0.0187 cbm
Weight of product	1.938 kg
Weight of case	1.94 kg
EAN of inner box	4038917141154
EAN case	4038917067140

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

