



Product description

Order number	74.5330.500
Product description	UD 3plus, Urine drainage system, 2 l, with holder, sterile, tube length: 100 cm, 1 piece(s)/blister
PZN number	16086848

Product characteristics

Colour of print/label	blue
Design	with holder
Application	Urine drainage
Volume of work	2 l

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**UD 3plus, Urine drainage system, 2 l, with holder,
sterile**

Page 2

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	1120
Depth of case	395 mm
Width of case	295 mm
Height of case	220 mm
Case volume	0.0256 cbm
Weight of product	0.191 kg
Weight of case	3.82 kg
EAN of inner box	6955824024333
EAN case	6955824024180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT