



Product description

Order number	74.5330.520
Product description	UD 3plus, Urine drainage system, 2 l, with holder, sterile, tube length: 120 cm, 1 piece(s)/blister
PZN number	16086831

Product characteristics

Colour of print/label	blue
Design	with holder
Application	Urine drainage
Volume of work	2 l

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	120 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	840
Depth of case	400 mm
Width of case	305 mm
Height of case	235 mm
Case volume	0.0287 cbm
Weight of product	0.179 kg
Weight of case	3.58 kg
EAN of inner box	6955824024340
EAN case	6955824024197

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

