



**Product Specification**  
**Urine drainage system, 2 l, with anti-reflux valve,  
with T drainage valve, sterile**

**Page 1**



**Product description**

Order number	74.5485.005
Product description	Urine drainage system, 2 l, with anti-reflux valve, with T drainage valve, sterile, tube length: 100 cm, 1 piece(s)/blister
PZN number	16356199

**Product characteristics**

Colour of print/label	blue
Design	with T drainage valve
Application	Urine drainage
Volume of work	2 l
Anti-reflux valve	yes

**Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine drainage system, 2 l, with anti-reflux valve,  
with T drainage valve, sterile**

**Page 2**

**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

**Purity & certification**

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2400
Depth of box	240 mm
Width of box	180 mm
Height of box	150 mm
Depth of case	395 mm
Width of case	295 mm
Height of case	290 mm
Case volume	0.0338 cbm
Weight of product	0.0844 kg
Weight of case	4.62 kg
EAN of inner box	6955824024418
EAN case	6955824024036

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**