



## Product description

|                     |   |
|---------------------|---|
| Order number        | 74.5684.001   |
| Product description | Urine measuring system exchange bag, 2 l, sterile, 1 piece(s)/blister |

## Product characteristics

|                       |                |
|-----------------------|----------------|
| Colour of print/label | blue           |
| Application           | Urine drainage |
| Volume of work        | 2 l            |

## Size

|               |     |
|---------------|-----|
| Sample volume | 2 l |
| Max. volume   | 2 l |

## Material & colours

|                   |                         |
|-------------------|-------------------------|
| Product material  | Polyvinylchloride (PVC) |
| Colour of product | transparent             |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

|                  |                     |
|------------------|---------------------|
| Product category | Medical device (MD) |
| CE certified     | CE supplier         |
| Purity standard  | sterile             |
| Sterilisation    | Gas sterilisation   |
| Batched          | yes                 |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 50            |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 50            |
| Piece(s) / pallet             | 4050          |
| Depth of box                  | 260 mm        |
| Width of box                  | 185 mm        |
| Height of box                 | 15 mm         |
| Depth of case                 | 395 mm        |
| Width of case                 | 195 mm        |
| Height of case                | 200 mm        |
| Case volume                   | 0.0154 cbm    |
| Weight of product             | 0.0451 kg     |
| Weight of case                | 3.02 kg       |
| EAN of inner box              | 6955824024432 |
| EAN case                      | 6955824024050 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

