



Product description

Order number	74.5684.002
Product description	Urine measuring system exchange bag, 2 l, with drainage valve, sterile, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	2 l

Size

Sample volume	2 l
Max. volume	2 l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine measuring system exchange bag, 2 l, with
drainage valve, sterile**

Page 2

Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	4050
Depth of case	395 mm
Width of case	245 mm
Height of case	200 mm
Case volume	0.0194 cbm
Weight of product	0.0401 kg
Weight of case	2.77 kg
EAN of inner box	6955824024449
EAN case	6955824024067

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP
support and is valid without signature.



SARSTEDT