



Product description

Order number	74.5767.001
Product description	UroStretch, size: S, material: PL
Aid number	15.25.11.1034
PZN number	8406948

Material & colours

Product material	polyester (PL)
Colour of product	white

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	10
Piece (s) / outer case	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
UroStretch, size: S

Page 2

Piece(s) / pallet	2340
Depth of box	283 mm
Width of box	173 mm
Height of box	135 mm
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0.0066 cbm
Weight of product	0.025 kg
Weight of case	0.37 kg
EAN of inner box	5703936100107
EAN case	5703936100145

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT