



**Product Specification**  
**Urine bag, 2 l, with anti-reflux valve, with drainage valve, sterile**

**Page 1**



**Product description**

Order number	74.5228.005
Product description	Urine bag, 2 l, with anti-reflux valve, with drainage valve, sterile, tube length: 90 cm, 1 piece(s)/blister
PZN number	16333726

**Product characteristics**

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	2 l
Anti-reflux valve	yes

**Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

**Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	100
Piece(s) / pallet	4800
Depth of case	305 mm
Width of case	390 mm
Height of case	255 mm
Case volume	0.0303 cbm
Weight of product	0.042 kg
Weight of case	4.5 kg
EAN of inner box	6971243661041
EAN case	6971243660808

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

