



## **Product description**

Order number	97.8720.500
Product description	SAHARA 4 100/240 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 4 storage bags, operation via touchscreen, automatic aggregate status recognition by heating zone and ability to switch blood products conveniently while warming device is running

## **Product characteristics**

Application	Warming
Power consumption during operation	1.000.00 W
Frequency	50 Hz - 60 Hz

## **Size**

Width of product	574 mm
Height of product	348 mm
Length of product	554 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	Medical device, CE
CE certified	CE
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.384 cbm
Weight of product	23.7 kg
Weight of case	27.7 kg
EAN of inner box	4038917968737
EAN case	4038917968744

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

