



Product description

Order number	79.8710.575
Product description	Paper roll, Consumables for protocol printer, white

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	7.5 cm
Length of product	8.3 cm

Material & colours

Product material	paper
Colour of product	white

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	10
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	50
Piece(s) / pallet	3000
Depth of case	400 mm
Width of case	400 mm
Height of case	170 mm
Case volume	0.0272 cbm
Weight of product	0.269 kg
Weight of case	14.15 kg
EAN of inner box	4038917975865
EAN case	4038917975872

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

