



Product description

| | |
|---------------------|--|
| Order number | 79.8710.576 |
| Product description | Ink ribbon, Consumables for protocol printer SP542MD, violet |

Product characteristics

| | |
|-------------|---------|
| Application | Warming |
|-------------|---------|

Size

| | |
|-------------------|-------|
| Width of product | 7 cm |
| Length of product | 12 cm |

Material & colours

| | |
|-------------------|--------|
| Colour of product | violet |
|-------------------|--------|

Purity & certification

| | |
|------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 144 |
| Piece(s) / pallet | 4608 |
| Depth of box | 120 mm |
| Width of box | 70 mm |
| Height of box | 30 mm |
| Depth of case | 480 mm |
| Width of case | 400 mm |
| Height of case | 230 mm |
| Case volume | 0.0442 cbm |
| Weight of product | 0.0311 kg |
| Weight of case | 7.5 kg |
| EAN of inner box | 4038917975902 |
| EAN case | 4038917975919 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

