



### **Product description**

Order number	97.8710.501
Product description	Collecting tray, for collecting plasma leaking from defective storage bags

### **Product characteristics**

Application	Warming
-------------	---------

### **Size**

Width of product	26 cm
Length of product	48 cm

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	72
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	1.5 kg
Weight of case	1.8 kg
EAN of inner box	4038917975827
EAN case	4038917975834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

