



Product description

Order number	97.8710.550
Product description	Module linfusion heater, for warming infusion solutions, tubes, instruments and contrast agent in the SAHARA-III

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	4.5 cm
Length of product	6 cm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of case	140 mm
Width of case	130 mm
Height of case	60 mm
Case volume	0.0011 cbm
Weight of product	0.03 kg
Weight of case	0.1 kg
EAN of inner box	4038917968379
EAN case	4038917968386

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

