



## **Product description**

Order number	97.8710.600
Product description	SAHARA-TSC 230 V, for thawing and warming blood products in plastic bags before transfusion, suitable for cryopreserved stem cells, incl. adaptation compress and warming tray

## **Product characteristics**

Application	Warming
Supply voltage	230 V
Power consumption during operation	655 W
Frequency	50 Hz - 60 Hz

## **Size**

Width of product	32 cm
Height of product	32.5 cm
Length of product	49.3 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	Medical device equipment, CE
CE certified	CE - manufacturer's self-declaration
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of case	800 mm
Width of case	460 mm
Height of case	470 mm
Case volume	0.173 cbm
Weight of product	16 kg
Weight of case	21 kg
EAN of inner box	4038917968539
EAN case	4038917968546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

