



### **Product description**

Order number	97.8710.620
Product description	Warming shell, Replacement part for storing stem cell preparations

### **Product characteristics**

Application	Warming
-------------	---------

#### **Size**

Width of product	18 cm
Length of product	27 cm

### **Material & colours**

Product material	aluminium
------------------	-----------

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0096 cbm
Weight of product	2.45 kg
Weight of case	2.95 kg
EAN of inner box	4028339720525
EAN case	4028339720518

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

