



Product description

Order number	79.8400.203
Product description	HAEM-o-TRANS set, collection and retransfusion of unwashed wound and drainage blood, set components: Y drain connector, suction bellows with anti-reflux barrier, HAEM-o-TRANS blood collection bag 700 ml

Product characteristics

Application	Autotransfusion
Connection	HAEM-o-TRANS

Purity & certification

Product category	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	Diethylhexylphthalate (DEHP)
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	20
Piece(s) / inner box	20
Piece(s) / case	20
Case volume	0,0562 cbm
Weight of product	0,165 kg
Weight of case	3,96 kg
EAN of inner box	4028339430226
EAN case	4028339430219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

