



Product description

Order number	79.8400.220
Product description	Blood collection bag, nominal volume: 700 ml, connection: HAEM-o-TRANS, pore size of pre-filter: 175 µm, with anti-reflux barrier, sterile

Product characteristics

Application	Autotransfusion
Filter	yes
Pore size	175
Connection	HAEM-o-TRANS

Size

Sample volume	700 ml
Width of product	19 cm
Length of product	25 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE0482 (Medcert)
Purity standard sterile
Sterilisation Gas sterilisation
Pyrogen-free yes
Plasticiser Diethylhexylphthalate (DEHP)
Batched yes

Packaging

Minimum order qty. 50
Piece(s) / inner box 10
Piece(s) / case 50
Piece(s) / pallet 3200
Case volume 0,0282 cbm
Weight of product 0,0466 kg
Weight of case 2,98 kg
EAN of inner box 4028339430325
EAN case 4028339430318

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

