



Product description

Order number	79.8400.230
Product description	HAEM-o-TRANS Transfusion device, with integrated cascade filter, 175 µm, 40 µm, 10 µm

Product characteristics

Application	Autotransfusion
-------------	-----------------

Purity & certification

Product category	Medical device
CE certified	CE 0482 (MedCert)
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	Diethylhexylphthalate (DEHP)
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	10
Type of smallest subpackaging	blister
Piece(s) / inner box	10
Piece (s) / outer case	80
Piece(s) / pallet	1600
Case volume	0.0616 cbm
Weight of product	0.0571 kg
Weight of case	5.33 kg
EAN of inner box	4028339430424
EAN case	4028339430417

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

