



### **Product description**

Order number	79.8400.270
Product description	Drainage bag, nominal volume: 700 ml, connection: HAEM-o-TRANS, with anti-reflux barrier, sterile

### **Product characteristics**

Application	Autotransfusion
Connection	HAEM-o-TRANS

### **Size**

Sample volume	700 ml
Width of product	20 cm
Length of product	23 cm

### **Material & colours**

Colour of product	transparent
-------------------	-------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category	Medical device
CE certified	CE0482 (Medcert)
Purity standard	sterile
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	yes
Batched	yes

### **Packaging**

Minimum order qty.	10
Piece(s) / inner box	10
Piece(s) / case	50
Case volume	0,0282 cbm
Weight of product	0,0258 kg
Weight of case	1,97 kg
EAN of inner box	4028339430523
EAN case	4028339430516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

