



Product description

Order number 97.8100.002
Product description Cover, Replacement part for K pressure infusion cuff

Product characteristics

Application Transfusion and infusion
Suitable for Pressure infusion cuff

Size

Max. volume 1000 ml
Width of product 14.5 cm
Height of product 2 cm
Length of product 27.5 cm

Material & colours

Colour of product red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Batched | no |

Packaging

| | |
|-------------------------|---------------|
| Minimum order qty. | 1 |
| Type of inner packaging | case |
| Piece(s) / inner box | 1 |
| Piece(s) / case | 5 |
| Case volume | 0.0039 cbm |
| Weight of product | 0.05 kg |
| Weight of case | 0.42 kg |
| EAN of inner box | 4028339310320 |
| EAN case | 4028339310313 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

