



Product description

Order number	97.8100.003
Product description	Manometer, Replacement part for the Classic Series pressure infusion cuff

Product characteristics

Application	Transfusion and infusion
-------------	--------------------------

Size

Width of product	4 cm
Height of product	2.5 cm
Length of product	7.5 cm

Material & colours

Colour of product	black
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 5
Case volume 0.0039 cbm
Weight of product 0.04 kg
Weight of case 0.37 kg
EAN of inner box 4028339310429
EAN case 4028339310412

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

