



### **Product description**

Order number	97.8100.004
Product description	Pressure ball, Replacement part for the Classic Series pressure infusion cuff, latex-free

### **Product characteristics**

Application	Transfusion and infusion
-------------	--------------------------

#### **Size**

Width of product	5 cm
Height of product	5 cm
Length of product	12 cm

### **Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Batched	no

### **Packaging**

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	5
Case volume	0.0039 cbm
Weight of product	0.07 kg
Weight of case	0.52 kg
EAN of inner box	4028339310528
EAN case	4028339310511

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

