



Product description

Order number	97.8100.010
Product description	Rubber bag, Replacement part for the Classic Series pressure infusion cuff, latex-free

Product characteristics

Application	Transfusion and infusion
-------------	--------------------------

Size

Width of product	12,5 cm
Height of product	1 cm
Length of product	25 cm

Material & colours

Product material	neoprene
Colour of product	black

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Latex-free	yes
Lot or serial number	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	5
Case volume	0,0039 cbm
Weight of product	0,11 kg
Weight of case	0,72 kg
EAN of inner box	4028339310627
EAN case	4028339310610

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

